

# Spectrums of activity and applications

	Disinfection	Cleaning	VAH listed	IHO virucide list	certified in acc. with EU standards	Bactericidal	Tuberculocidal	Yeasticidal	Fungicidal	Limited virucidal	Limited virucidal plus	Virucidal	Concentration / Reaction time in min.	No. of days product works for
<b>Instruments</b>														
ID 215 Enzymatic instrument cleaner	■	■			■	■		■		■			1%/2%/30/15	
<b>Surfaces</b>														
FD 322 Quick-acting disinfection	■		■	■	■	■	■	■	■	■	■		100%/1*	
FD 360 Cleaning and care of vinyl upholstery		■											ready to use /1*	
FD 366 sensitive wipes for disinfection of sensitive surfaces	■	■	■	■	■	■	■	■		■			ready to use /1*	
FD 366 sensitive top wipes for disinfection of sensitive surfaces	■	■	■	■	■	■	■	■		■			ready to use /1*	
<b>Special areas</b>														
Orotol® plus Disinfection of suction systems	■	■	■	■	■	■		■		■			2%/5*	
MD 555 cleaner														
Special detergent for suction units		■											5%/30-120	

\* bactericidal/yeasticidal.

For more information, please refer to the product information or visit [www.duerdental.com](http://www.duerdental.com)



DÜRR DENTAL OCEANIA PTY LTD.

7/42 Sabre Drive

Port Melbourne, Victoria 3207

Australia

E: [mina.farag@duerrdental.com](mailto:mina.farag@duerrdental.com)

T: +61 437 799 323

T: 1300 52 53 51

**DÜRR  
DENTAL**  
THE BEST, BY DESIGN

# System Hygiene from Dürr Dental



## ■ Orotol® plus Disinfection of suction systems



Foam-free liquid concentrate for the simultaneous disinfection, cleaning, deodorisation and care of all suction systems and amalgam separators. Good for all suction systems. Recommended by DentsplySirona, Planmeca, Morita, Belmont & Heka.



## ■ MD 555 cleaner Special cleaner for suction systems



Foam-free cleaning concentrate for all kinds of suction systems and amalgam separators. Recommended by DentsplySirona, Planmeca & Heka.



## ■ OroCup Care system

User-friendly, closed care system for the disinfection and cleaning of suction systems and amalgam separators. Use with Orotol® plus & MD 555 cleaner.



## ■ ID 215 Enzymatic Cleaner



Concentrate for manual, enzymatic instrument cleaning prior to disinfection.



## ■ Hygobox

Sturdy and robust transport and disinfection box with a 3 litre capacity



## ■ FD 322 Quick-acting disinfection



Fast-acting, ready-to-use solution for the disinfection of surfaces. Recommended by Planmeca, KAVO & Heka.



## ■ FD 366 Sensitive wipes and top wipes Disinfection of sensitive surfaces



Ready-to-use disinfection and cleaning wipes for sensitive surfaces. Recommended by A-dec, DentsplySirona, Morita, Belmont & Heka.

**Available as:** Dispenser box with 100 wipes  
**Wipe size:** 20 x 30 cm



## ■ FD 366 Sensitive disinfection of sensitive surfaces



Alcohol-reduced, ready-to-use solution for the disinfection and cleaning of especially sensitive surfaces. Recommended by A-dec, DentsplySirona, Morita, Belmont & Heka.



## ■ FD 360



Cleaning and care of vinyl upholstery  
Recommended by Planmeca, DentsplySirona, Morita & Heka.



## ■ Hygoclave 40/50 Series

The Hygoclave 40 & 50 Series are a premium autoclave made by Dürr Dental. It includes a high-resolution 5" touch screen, flexible water supply options, and cloud connectivity with the 50 Series. Recommended by A-dec.