



# System Hygiene from Dürr Dental



## ■ ID 215 Enzymatic Cleaner

Concentrate for manual, enzymatic instrument cleaning prior to disinfection.



## ■ FD 322 Quick-acting disinfection

Fast-acting, ready-to-use solution for the disinfection and cleaning of surfaces. Recommended by Planmeca, KAVO & Heka.



## ■ FD 333 Quick-acting disinfection

Ready-to-use solution for the quick-acting disinfection of the surfaces of medical devices and medical inventory. Recommended by Planmeca, Morita, Belmont & Heka.



## ■ FD 333 top wipes Quick-acting disinfection

Alcohol-based, ready-to-use disinfection and cleaning wipes for smaller surfaces of medical devices. Recommended by Belmont, Morita, Planmeca & Heka.



## ■ FD 360

Cleaning and care of vinyl upholstery. Recommended by Planmeca, DentsplySirona, Morita & Heka.



## ■ FD 366 Sensitive disinfection of sensitive surfaces

Alcohol-reduced, ready-to-use solution for the disinfection and cleaning of especially sensitive surfaces. Recommended by A-dec, DentsplySirona, Morita, Belmont & Heka.



## ■ FD 366 Sensitive wipes and top wipes Disinfection of sensitive surfaces

Ready-to-use disinfection and cleaning wipes for sensitive surfaces. Recommended by A-dec & Dentsply Sirona.

Available as: Dispenser box with 100 wipes  
Wipe size: 20 x 30 cm



## ■ MD 550 Spittoon bowl cleaner

Foam-free, ready-to-use solution for intensive, residue-free cleaning of spittoons.

Recommended by Dentsply Sirona, Morita & Heka.

## ■ Orotol<sup>fl</sup> plus Disinfection of suction systems



Foam-free liquid concentrate for the simultaneous disinfection, cleaning, deodorisation and care of all suction systems and amalgam separators. Good for all suction systems.

Recommended by DentsplySirona, Planmeca, Morita, Belmont & Heka.



## ■ MD 555 cleaner Special cleaner for suction systems

Foam-free cleaning concentrate for all kinds of suction systems and amalgam separators.

Recommended by DentsplySirona, Planmeca & Heka.

## ■ OroCup Care system



User-friendly, closed care system for the disinfection and cleaning of suction systems and amalgam separators. Use with Orotol Plus & MD555.